

# Anspach XMax. Pneumatic high speed performance and power in a variety of applications.

Smooth performance

Continuous lubrication system

Effortless assembly



High Speed  
Power Tools

**Anspach XMax.** Pneumatic high speed performance and power in a variety of applications.


## Range of applications



ENT



Neurosurgery



Spine surgery

## High performance handpiece with an advanced design

- Ergonomic design
- Lightweight handpiece (75 g), optimal for long procedures
- Powerful system (90 watts), ideal for demanding applications
- Very quiet system
- Very precise thanks to a smooth performance with minimal start-up kick
- Effortless assembly thanks to a quick connect/disconnect coupling



High flexibility and balance thanks to an angled swivel elbow: the motor moves independently from the hose

Innovative coupling system for quick and safe attachment exchange with no tool needed

## Foot control for continuously variable and precise speed control

### Foot control

- Continuous lubrication integrated in the pedal for an increased reliability of the system
- Easy set-up



## Comprehensive line of attachments and accessories for maximum versatility

- Both general attachments for standard applications and surgery specific attachments available
- Quick connect/disconnect system for a hassle free attachment and burr changing during surgery
- Color-coding on attachment and burr for more simplicity
- Identical attachments as for the eMax 2 Plus

### Straight attachments

Straight



Straight Heavy Duty (HD)



Craniotomes



### Compact Speed Reducer



### Angled attachments

Quick Disconnect Angle Attachments (QD)



Minimal Access Attachments (MA)



**Anspach XMax.** Pneumatic high speed performance and power in a variety of applications.

## Ordering information

This brochure does not present all the existing Anspach articles. For a complete overview about the XMax system, please refer to the Synthes e-catalogue.

### Handpiece

XMAX	XMax Handpiece
------	----------------

### Foot controls

AUTOLUBE-III	Pneumatic Foot Control
AUTOLUBE-III-IRR	Pneumatic Foot Control for Irrigation System
AUTOLUBE-III-MA7	Pneumatic Foot Control with MA7

### Attachments

SHORT	5 cm Short Attachment. Use with XMax, microMax, and eMax Systems.
MEDIUM	8 cm Medium Attachment. Use with XMax, microMax, and eMax Systems.
LONG	11 cm Long Attachment. Use with XMax, microMax, and eMax Systems.
MIA16	16 cm Minimally Invasive Attachment. Use with XMax, microMax, and eMax Systems.
SHORT-HD	5 cm Short Heavy Duty Attachment. Use with XMax, microMax, and eMax Systems.
MEDIUM-HD	8 cm Medium Heavy Duty Attachment. Use with XMax, microMax, and eMax Systems.
LONG-HD	12.4 cm Heavy Duty Long Attachment. Use with XMax, microMax, and eMax Systems.
XL-HD	20.2 cm Heavy Duty Extra Long Attachment. Use with XMax, microMax, and eMax Systems.
ORANGE-45	45 Contra Angle Attachment, Multiple Angle. Use with XMax, microMax, and eMax Systems.
ORANGE-90	90 Right Angle Attachment. Use with XMax, microMax, and eMax Systems.
CRANI-A	6.5 cm Adult Craniotome Attachment. Use with XMax, microMax, and eMax Systems.
CRANI-A-01	6.5 cm Adult Craniotome, Thin Foot Plate Attachment. Use with XMax, microMax, and eMax Systems.
CRANI-A-R	6.5 cm Adult Rotating Craniotome Attachment. Use with XMax, microMax, and eMax Systems.
CRANI-L	7.5 cm Large Craniotome Attachment. Use with XMax, microMax, and eMax Systems.
CRANI-L-R	7.5 cm Large Rotating Craniotome Attachment. Use with XMax, microMax, and eMax Systems.
CRANI-P	6.5 cm Pediatric Craniotome Attachment. Use with XMax, microMax, and eMax Systems.

ADG	7.2 cm Adjustable Drill Guide. Adjusts 3–6 mm increments. Use with XMax, microMax, and eMax Systems.	
CDA	6.8 cm Controlled Depth Attachment. Adjusts 1–4 mm increments. Use with XMax, microMax, and eMax Systems.	
MA-DRIVER	Minimal Access Straight Driver. Use with XMax, microMax, and eMax Systems.	
MA-D20	Minimal Access 20 Angle Driver. Use with XMax, microMax, and eMax Systems.	
Bearing sleeve	Length (cm)	Type
MA-10S	10	straight
MA-10C	10	curved
MA-15S	15	straight
MA-15ST	15	straight tapered
MA-15C	15	curved
MA-19ST	19	straight tapered
QD8	8 cm Angle Attachment. Use with XMax, microMax, and eMax Systems.	
QD11	11 cm Angle Attachment. Use with XMax, microMax, and eMax Systems.	
QD14	14 cm Angle Attachment. Use with XMax, microMax, and eMax Systems.	

### Other attachments and couplers

CSR60	Perforator Driver with Hudson End, 60:1 Ratio. Use with XMax, microMax, and eMax Systems.
SA-JACOBS	Small Attachment, Keyless Jacobs Chuck. Use with all Anspach Systems.
SA-JLATCH	Small Attachment, J-Latch. Use with all Anspach Systems.

### Other

ACB	Attachment Cleaning Brushes (One 7" Brush and One 12" Brush)
B-OIL	Instrument Oil and One Replacement Filter for AUTOLUBE-I
FILTER-II	Five Replacement Filters for AUTOLUBE-II or AUTOLUBE-III

### Irrigation

For ordering information about the Irrigation Clips, Tubings or Pump, please refer to the Synthes e-catalogue.

### Cutting Tools

For ordering information about the XMax Cutting Tools, please refer to the Cutting Burr Charts (035.000.500, 035.000.505).

### ECREP

Anspach Europe Ltd.  
Aston Court, Kingsmead Business Park  
Frederick Place, High Wycombe  
HP11 1LA, United Kingdom



The Anspach Effort, Inc.  
4500 Riverside Drive  
Palm Beach Gardens, FL 33410  
USA

Distributor:  
Synthes GmbH  
Eimattstrasse 3  
CH-4436 Oberdorf  
www.synthes.com



All technique guides are available as PDF files at  
www.synthes.com/lit

